

SP SUPERTEST 450 - MIX TESTOSTERONE (TESTOSTERONE 5) SP LABS DESCRIPTION

With great joy we present you a super new product from the company **SP Laboratories - Supertest 450** .

What is Supertest 450?

This is as much as 5 different testosterone esters with different speeds on and half-life in one product! 450 mg. of active substance per 1 ml.

- Test acetate 32 mg.
- Test propionate 73 mg.
- Test phenylpropionate 73 mg.
- Test cypionate 125 mg.
- Test decanoate 147 mg.

Supertest 450 can be loosely compared to Sustanon. The similarity is expressed in the duration of exposure, and in the presence of two identical esters. The main difference is that **Supertest 450** not inherent presence of "hormonal pits", and this is its main advantage lies, and the optional ether acetate, cypionate.

Why you should **buy 450 SP Supertest**

- Turning to work within a few minutes after injection
- Smooth switching from short to long ester keeps stable "working" the concentration of testosterone in the body.
- Savings of 30-40%. In 10 mL vial, in fact, equivalent to the contents of two bottles of any other testosterone. Enough time.
- Super quality of JV Lab
- A powerful effect on the anabolic processes in the body

Supertest 450 application, course reviews

Supertest 450 on the duration of exposure is the "long-playing". 1 After injection, the drug will work for almost a month. And based on this, short term courses are not recommended before 10 weeks.

Course on superteste 450 should last 10 weeks, for what would be the fullest to feel the effect of the drug.

Like any testosterone, Supertest 450 can be combined with any drugs, according to the scheme + anabolic androgen. Thus, you will achieve a synergistic effect of the drug and greatly increase the effectiveness of the course.

Supertest side effects post-cycle therapy

Supertest susceptible to flavor, just like any other testosterone, so if you plan on Superteste course, it is necessary to 2 weeks of the course start receiving anastrozole, leveling the side effects associated with the conversion to estrogen SP Supertest.

Also be sure to use gonadotropin for 3 weeks before the end of the course. The use of gonadotropin should begin with the last injection supertesta.

Preparations on the FCT are selected depending on with which the drug was combined supertest. If the course is not used nandrolone trenbolone, then as a primary prepa FCT fit tamoxifen . If you used the above mentioned products, it is necessary to conduct pips Clomid. PBC should begin after the complete excretion of the drug - 3.5 weeks.

By **upit and order Supertest 450 by SP Laboratories** you can at the lowest price in the USA. We guarantee the original drug.